


Red tech official for radiographer

10 MCQS
X-RAY TUBE CONSTRUCTION
& WORKING

1. What is the typical material used for the filament in an X-ray tube?

- A. Copper
- B. Molybdenum
- C. Tungsten
- D. Lead




Correct answer: C. Tungsten

2. Which part of the X-ray tube rotates to dissipate heat?

- A. Filament
- B. Focusing cup
- C. Anode disk
- D. Cathode cup

Correct answer: : C. Anode disk



3. The angle of the anode target in most general-purpose X-ray tubes is:

A. 5°

B. $10-20^{\circ}$

C. $30-40^{\circ}$

D. $45-60^{\circ}$

Correct answer: B. $10-20^{\circ}$

4. The line focus principle is used in X-ray tubes to:
- A. Increase anode rotation speed
 - B. Produce a smaller effective focal spot
 - C. Improve beam filtration
 - D. Prevent tube arcing

Correct answer: B. Produce a smaller effective focal spot

5. The envelope of the X-ray tube is usually made of:

A. Plastic

B. Steel

C. Lead

D. Pyrex glass or metal

Correct answer: D. Pyrex glass or metal

6. Which component allows rotation of the anode?

- A. Rotor and stator assembly
- B. Filament transformer
- C. Grid circuit
- D. Cathode bias unit

Correct answer: A. Rotor and stator assembly

7. The term “space charge effect” refers to:

A. Anode overheating

B. X-ray beam hardening

C. Electron cloud repelling further emission

D. Incomplete beam collimation

Correct answer: C. Electron cloud repelling further emission

8. The primary cause of X-ray tube failure is:

A. Grid failure

B. Filament evaporation

C. Window cracking

D. Focusing cup distortion

Correct answer: B. Filament evaporation

9. The purpose of the tube housing is:

A. Focus the electron beam

B. Generate the X-rays

C. Limit leakage radiation and support cooling

D. Vacuum regulation

Correct answer: C. Limit leakage radiation and support cooling

10. Dual-focus X-ray tubes are designed to:

- A. Use two rotating anodes
- B. Provide large and small focal spots
- C. Operate at higher kVp
- D. Allow dual-energy CT

Correct answer: B. Provide large and small focal spots